

Callis Ortho

For motion analysis and orthopedic training



The Callis Ortho has been specifically developed for motion/ gait analysis and orthopedic training where the accuracy and durability of our slat-belt running surface provide an unparalleled gait experience.

Complete with reverse running mode, open and shortened hand rails, downhill mode and permanent white marking in the middle of the belts, the Callis Ortho fits with every orthopedic and motion lab.

Technical data

Running surface:	Low-impact, shock-absorbing slat-belts
Running surface dimensions:	61x20 in / 155x50 cm
Weight capacity:	< 353 lbs / 160 kg
Speed range:	0-13.7 mph / 0-22 km/h in 0.1 mph increments
Dimensions (LxWxH):	68x30x51 in / 172x76x130 cm (without add-ons)
Weight:	ca. 353 lbs / 160 kg (without add-ons)
Noise level:	<75 db (A)
Motor:	1.3 kW 230 V AC, 50/60 Hz, 6.3 A

All specifications are subject to change without notice.

[Type here]

Benefits for motion analysis

Reverse running direction
Downhill running mode
Permanent white marking in the middle of the surface

Suitable for every requirement

0-20% incline range
Max. weight capacity increase up to 250kg
Removable handrail keypad
Additional emergency button on handrail
Serial interface RS232, galvanic isolation



UK Distributors – Mar Systems Ltd

Unit E6, Ascot Business Park
Lyndhurst Rd, Ascot, Berkshire
SL5 9HS

Tel 01344 623 883

[Type here]

www.mar-systems.co.uk

info@mar-systems.co.uk